



# Recombinant Mouse Proto-oncogene c-Rel (Rel)

<b>Product Code</b>	CSB-EP019552MO-B
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P15307
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	ASSGYNPYV EIIEQPRQRG MRFYKCEGR SAGSIPGERS TDNNRTYPSV QIMNYYGK GK IRITLVTKND PYKPHPHDLV GKDCRDGYYE AEFGPERRPL FFQNLGIRCV KKKEVKGAI LRISAGINPF NVGEQQLLDI EDCDLNVVRC VFMFFLPDED GNFTTALPPI VSNPIYDNRA PNTAELRICR VNKNCGSVRG GDEIFLLCDK VQKDDIEVRF VLNDWEARGV FSQADVHRQV AIVFKTPPYC KAILEPVTVK MQLRRPSDQE VSESMDFRYL PDEKDAYGNK SKKQKTTLIF QKLLQDCGHF TEKPRTAPLG STGEGRFIKK ESNLFSHGTV LPEMPRSSGV PGQAEPYYSS CGSISGLPH HPPAIPVAH QPTSWPPVTH PTSHPVSTNT LSTFSAGTLS SNSQGILPFL EPGVSDLSA SNSCLYNPDD LARMETPSMS PTDLYISDV NMLSTRPLSV MAPSTDGMGD TDNPRLV SIN LENPSCNARL GPRDLRQLHQ MSPASLSAGT SSSSVFVSQS DAFDRSNFSC VDNGLMNEPG LSDDANNPTF VQSSHYSVNT LQSEQLSDPF TYGFFKI
<b>Source</b>	E.coli
<b>Target Names</b>	Rel
<b>Protein Names</b>	Recommended name: Proto-oncogene c-Rel
<b>Expression Region</b>	2-587
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.