



# Recombinant Arabidopsis thaliana GDSL esterase/lipase At1g54000 (At1g54000)

<b>Product Code</b>	CSB-BP625342DOA
<b>Abbreviation</b>	At1g54000
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q1H583
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	Q NSPVVALFTF GDSNFDAGNK QLTCLKLVAQ GFWPYGKSRD DPNGKFS DGL ITPDFLAKFM KIPLAIAPAL QPNVNVSRGA SFAVEGATLL GAPVESMTLN QQVKKFNQMK AANWNDDFVA KSVFMIIYIGA NDYLNFTKNN PTADASAQQA FVTSVTNKLK NDISALYSSG ASKFVIQTLA PLGCLPIVRQ EYNTGMDQCY EKLNDLAKQH NEKIGPMLNE MARNSPASAP FQFTVDFDYN AVLTRTQRNQ NFRFFVTNAS CCGVGSHDAY GCGLPNVHSK LCEYQRSFLF FDGRHNSEKA QEMFAHLLFG ADTNVVQPMN VRELTVYPVD EPMREFWVPP TPATVHASDS SSSTSRGYEY Y
<b>Source</b>	Baculovirus
<b>Target Names</b>	GLL22
<b>Protein Names</b>	Recommended name: GDSL esterase/lipase At1g54000 EC= 3.1.1.- Alternative name(s): Extracellular lipase At1g54000
<b>Expression Region</b>	30-391
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.