



Recombinant Arabidopsis thaliana GDSL esterase/lipase At1g54000 (At1g54000)

Product Code	CSB-YP625342DOA
Abbreviation	At1g54000
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q1H583
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	Q NSPVVALFTF GDSNFDAGNK QTLTKTLVAQ GFWPYGKSRD DPNGKFS DGL ITPDFLAKFM KIPLAIAPAL QPNVNVSRGA SFAVEGATLL GAPVESMTLN QQVKKFNQMK AANWNDDFVA KSVFMIYIGA NDYLNFTKNN PTADASAQQA FVTSVTNKLK NDISALYSSG ASKFVIQTLA PLGCLPIVRQ EYNTGMDQCY EKLNDLAKQH NEKIGPMLNE MARNSPASAP FQFTVDFDYN AVLTRTRQRNQ NFRFFVTNAS CCGVGSHDAY GCGLPNVHSK LCEYQRSFLF FDGRHNSEKA QEMFAHLLFG ADTNVVQPMN VRELTVYPVD EPMREFWVPP TPATVHASDS SSSTSRGYEY Y
Source	Yeast
Target Names	GLL22
Protein Names	Recommended name: GDSL esterase/lipase At1g54000 EC= 3.1.1.- Alternative name(s): Extracellular lipase At1g54000
Expression Region	30-391
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.