



Recombinant Danio rerio UPF0471 protein C1orf63 homolog (zgc:136474)

Product Code	CSB-BP626477DIL
Abbreviation	zgc:136474
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q1ECZ9
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MKTEASPGRL HEDVKLIFDK KAPSGRSSSC SSSSSSSSGS SYSSSRGSR SRRSRTSSSS RSSSSGSSSS SNSRSRSPR CSGRVHCRHR HRSPRRYRA RRSYSPSPE SSSRRRHYYR RSRSPYRYSR RYRRSPSRSR SRSPPHWRGS RFIGRYRCRF SQSPRRSPRP YRSRSRER SIRLSLEEK YLLNVAKANA ARILGVQNL LPESLKEMEQ QEERKRSSS DEEERVVDL AAQKTPAQVN GGAADDEAET AQTSPKRKQI VFSNNNAIAK PSSSPTLSDS KVTSRADIVP YKRPFQWVP IKGSSSANKR
Source	Baculovirus
Target Names	rsrp1
Protein Names	Recommended name: UPF0471 protein C1orf63 homolog
Expression Region	1-320
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.