



Recombinant Danio rerio Homeobox protein Hox-B13a (hoxb13a)

Product Code	CSB-MP627339DIL
Abbreviation	hoxb13a
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q1ECY2
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MTTSLVLNPR WADTVMFVYE>NNLDELKNME GLVGSGNFAA NQCRNLMAHS ALSGHPSSLV HGSSYPTVDV STSSSAESGK QCTPCPTVPQ ASSTGPIPYG YFGNSYYPCR MGRGSLKSCT QPSALSYTAE KYMDTPVTSE EYPTRAKEFA FYHSYSPYQ SMASYLDVSV VQTLGTGEPR HDSLLPMSY QPWALANGWG SQMYCSKDQG QAGHLWKSAL ADVVAHQHDG GSFRRGRKKR IPYTKVQLKE LEKEYAANKF ITKDKRRKIS AVTNLSERQI TIWFQNRVRK EKKFIAKVKN TAP
Source	Mammalian cell
Target Names	hoxb13a
Protein Names	Recommended name: Homeobox protein Hox-B13a
Expression Region	1-303
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.