



Recombinant Aedes aegypti Probable dynactin subunit 2 (Dmn)

Product Code	CSB-BP627929AXQ
Abbreviation	Dmn
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q1HQF2
Product Type	Recombinant Protein
Immunogen Species	Aedes aegypti (Yellowfever mosquito) (Culex aegypti)
Purity	>85% (SDS-PAGE)
Sequence	MADPKFQFLP YIAHDQPDVY ETPDAPESET SDFYDEEPAN ESIERLHIST KDSYNKFKGK YLTGNVDFSD RIGKKIRNGY DAVSGEWDLA GEGETETPIQ KCLRLQCEMN ELMEEITASQ ADATKTKEEK ASYEAVFDVV GTAKKVLESL KLEQAIGSET VAMGAETEA KLLTKIEEYK KSVGADPAKT ATELVYSSRI AELEHRLHEL EVAVGAKPEK ISRLAGSAGA GNLEAVQSI SAKAALLQPQ QLDLIETRLG NLLLKMAAID EKSAATGQDA NREQKILELY EIAKTTEPIV QILPDMLNRM QTLESLHKYA TNFSKLFAEL ETTQASILKG IASNKTLTGT VQEAFANLE NVNKEVKKLD ERMKTLQEKV K
Source	Baculovirus
Target Names	DCTN2-p50
Protein Names	Recommended name: Probable dynactin subunit 2
Expression Region	1-381
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.