



Recombinant Takifugu rubripes Src kinase-associated phosphoprotein 1 (skap1)

Product Code	CSB-EP629278TKY
Abbreviation	skap1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q1KKW7
Product Type	Recombinant Protein
Immunogen Species	Takifugu rubripes (Japanese pufferfish) (Fugu rubripes)
Purity	≥85% (SDS-PAGE)
Sequence	MEAVDLLANV SVTGLLLFWG DCEFFVSEIL YDENLSENAQ ETRKVLLNNF RVVHVRNPQE FPFPSDYKEE DGSDDNRSSS LGRSAQSDDA SLASDNQDDG IPEYFGDIPP VAAQDIVNVL KQGYLEKKRK DHSFFGSEWQ KRWCVLNNLV FYYFGSDKDK QQKGSFYISD YSVQLVTNLR KDSRKTSCFE FFAPGRRPFQ FTAGSPQEAKEWVDQIKIVL RDLSSSTVIPV DDEEEEEEEEE ETYDDIEGEG GPPLPQPLSG TWGRGGDTGA ADEEDEDIYE VLPEESPDSA DGSMERNKPK EYANYQQLW DCSADEPDEL PFQRGDLIYI ISKEYNIYGW WVGELNGAVG IVPKDFLHPA YIL
Source	E.coli
Target Names	skap1
Protein Names	Recommended name: Src kinase-associated phosphoprotein 1 Alternative name(s): Src family-associated phosphoprotein 1
Expression Region	1-363
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.