



# Recombinant Human Heterogeneous nuclear ribonucleoprotein U-like protein 2 (HNRNPUL2)

<b>Product Code</b>	CSB-MP629282HU
<b>Abbreviation</b>	HNRNPUL2
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q1KMD3
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MEVKRLKVTE LRSELQRRGL DSRGLKVDLA QRLQEALDAE MLEDEAGGGG AGPGGACKAE PRPVAASGGG PGGDEEEDEE EEEDEEALL EDEDEEPPPA QALGQAAQPP PEPPEAAAME AAAEPDASEK PAEATAGSSG VNGGEEQGLG KREDEPEER SGDETPGSEV PGDKAAEEQG DDQDSEKSKP AGSDGERRGV KRQRDEKDEH GRAYYEFREE AYHSRSKSPL PPEEEAKDEE EDQTLVNLDT YTSDLHFQVS KDRYGGQPLF SEKFPTLWSG ARSTYGVTKG KVCFEAKVTQ NLPMKEGCTE VSLLRVGWSV DFSRPQLGED EFSYGFDRG LKAENGQFEE FGQTFGENDV IGCANFETE EVELSFSKNG EDLGVAFWIS KDSLADRALL PHVLCKNCVV ELNFGQKEEP FFPPPEEFVF IHAVPVEERV RTAVPPKTIE ECEVILMVGL PGSGKTQWAL KYAKENPEKR YNVLGAETVL NQMRMKGLEE PEMDPKSRDL LVQQASQCLS KLVQIASRTK RNFILDQCNV YNSGQRRKLL LFKTFSRKVV VVVPNEEDWK KRLELRKEVE GDDVPESIML EMKANFSLPE KCDYMDEVTY GELEKEEAQP IVTKYKEEAR KLLPPSEKRT NRRNNRNKRN RQNRSRGQGY VGGQRRGYDN RAYGQQYWGQ PGNRGGYRNF YDRYRGDYDR FYGRDYEYNR YRDYYRQYNR DWQSYYYHHP QDRDRYRNY YGYQGYR
<b>Source</b>	Mammalian cell
<b>Target Names</b>	HNRNPUL2
<b>Protein Names</b>	Recommended name: Heterogeneous nuclear ribonucleoprotein U-like protein 2 Alternative name(s): Scaffold-attachment factor A2 Short name= SAF-A2
<b>Expression Region</b>	1-747
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.