



Recombinant Danio rerio UPF0583 protein C15orf59 homolog (si:dkey-105e17.1)

Product Code	CSB-MP629723DIL
Abbreviation	si:dkey-105e17.1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q1L899
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MSQSRAPARE PSETPSQREQ IRSHMKMVIQ QLEGILKELK DVAHELREVV GQIDKLTSDL ELDLDADDWT VATASSTSSS ERGLCEAFRL DFLGQDSLSD SWDFCSFLES SSRRSARDDT KPPPTTASVY SQMNGGLPVP NGPLIITPDS SSEEASSSTH SQSQKTSRTA GTRERVRFSK KILYHALCCD DDEDEDEDED GRDEEEDKLD TDSERSPLAG SPVPPLPELY NCREPPAGGS HTIPRKDALN PGCRKLLRN SSTQTVSDKS TQTLIPYVPA KHRNKDL
Source	Mammalian cell
Target Names	si:dkey-105e17.1
Protein Names	Recommended name: UPF0583 protein C15orf59 homolog
Expression Region	1-287
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.