



Recombinant Schizosaccharomyces pombe Putative histone deacetylase complex subunit cti6 (cti6)

Product Code	CSB-MP629890SXV
Abbreviation	cti6
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q1MTR4
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MPSNSTPANE EKEEESKIQ KSPVDTPHTP NGSVSDNEEN ETSSTGEVTR CVCGIVESDD EASDGGLYIQ CDQCSVWQH G NCVGFADESE VPEVYYCEIC HPEFHKVYQR GRGAKQSKYL GNGKPIEASQ TEESSTPPS PATKKSSKQR LTMNSRDAAL DYEEYLAIK EKSLIPRRSR GRTSSKSLSP PAPQDESQGT EINLKQKIEE ENDEILEDK ESKDENEENK ETSTTNVAET DAPEEETVDT VEEIADEEKH SVKEESGEAS PQSSQQSTIT SISTTTRSTR KAKREAAAED KADLPAAVAP KPSKTRKVGG RRGKSSSNDN HRIPQLHPDG TFVETITKPK GLHSRITMTE MRRRVASMLE YIGHIQVEMA AQSAGNQSST KSSKEGPEEE KETLRMVDNL TRDLLHWEQR FSRT
Source	Mammalian cell
Target Names	cti6
Protein Names	Recommended name: Putative histone deacetylase complex subunit cti6
Expression Region	1-424
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.