



# Recombinant Arabidopsis thaliana Pentatricopeptide repeat-containing protein At1g28690, mitochondrial (PCMP-E34), partial

<b>Product Code</b>	CSB-MP629900DOA
<b>Abbreviation</b>	PCMP-E34
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q1PFQ9
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	S TFPLKQNVSS LSPAKYIAGA LQEHINSPAP KAGKKIHADI IKTFGFQPDLN ISIKLLILHL KCGCLSYARQ VFDELPKPTL SAYNYMISGY LKHGLVKELL LLVQRMSYSG EKADGYTLSM VLKASNSRGS TMILPRSLCR LVHARIKCD VELDDVLITA LVDTYVKSGK LESARTVFET MKDENVCCT SMISGYMNQG FVEDAEEIFN TTKVKDIVVY NAMVEGFSRS GETAKRSVDM YISMQRAGFH PNISTFASVI GACSVLTSHE VGQQVHAQIM KSGVYTHIKM GSSLLDMYAK CGGINDARRV FDQMKEKNVF SWTSMIDGYG KNGNPEEAL LFTRMKEFRI EPNYVTFLGA LSACSHSGLV DKGYEIFESM QRDYSMKPKM EHYACIVDLM GRAGDLNKAF EFARAMPERP DSDIWAALLS SCNLHGNVEL ASIAASELFK LNADKRPGAY LALSNVYASN DKWDNVSKIR EVMKRRRISK TIGRSWTSSE
<b>Source</b>	Mammalian cell
<b>Target Names</b>	PCMP-E34
<b>Protein Names</b>	Recommended name: Pentatricopeptide repeat-containing protein At1g28690, mitochondrial
<b>Expression Region</b>	20-520
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Partial
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at  $-20^{\circ}\text{C}/-80^{\circ}\text{C}$ . The shelf life of lyophilized form is 12 months at  $-20^{\circ}\text{C}/-80^{\circ}\text{C}$ .