



Recombinant Nitrobacter hamburgensis Fructose-1,6-bisphosphatase class 1 2 (fbp2)

Product Code	CSB-EP629943NAAF
Abbreviation	fbp2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q1QGD7
Product Type	Recombinant Protein
Immunogen Species	Nitrobacter hamburgensis (strain DSM 10229 / NCIMB 13809 / X14)
Purity	≥85% (SDS-PAGE)
Sequence	MSHELTLQQR LEVQAGPDPL HRVAVKEAVAA LAQAAIDISD LTCRGALAGI TGEAQGRNTD GDIQKDLVDV ADQIIRDALS ALPIAVLASE EMAELDILNP GAPISVAFDP LDGSSNINTN ISVGTIFSIM PTPPDASAAF TQPGSAQLAA GFVVYGPQTS LILTLGQGVD IFTLDRVERV FKLTGSMMQI PADATEFAIN TSNRRHWDLP VRAYIDECLM GADGPSGKNF NMRWIGSLVA EAFRILVRGG IFLYPGDARD GYEDGRLRLV YEAHPMAFII EQAGGGASTG RKRILDIVPG SLHQRVPLIM GSIKNVRRLE DMHTGPNVAL EANAPLFGHR GLFRV
Source	E.coli
Target Names	fbp2
Protein Names	Recommended name: Fructose-1,6-bisphosphatase class 1 2 Short name=FBPase class 1 2 EC= 3.1.3.11 Alternative name(s): D-fructose-1,6-bisphosphate 1-phosphohydrolase class 1 2
Expression Region	1-345
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.