



# Recombinant Escherichia coli Glycogen debranching enzyme (glgX)

<b>Product Code</b>	CSB-EP630059EGW-B
<b>Abbreviation</b>	glgX
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q1R5J5
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli (strain UTI89 / UPEC)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MTQLAIGKPT PLGAHYDGQG VNFTLFSAHA ERVELCVFDA NGQEHRYDLP GHSGDIWHGY LPDARPGLRY GYRVHGPWQP AEGHRFNPAK LLIDPCARQI DGEFKDNPLL HAGHNEPDYR DNASIAPKCV VVDHYDWED DAPPRMPWGC TIIYEAHVKG LTYLHPEIPV EIRGTYKALG HPVMINYLKQ LGITALELLP VAQFASEPRL QRMGLSNYWG YNPVAMFALH PAYACSPETA LDEFRDAIKA LHKAGIEVIL DIVLNHSAEL DLDGPLFLSLR GIDNRSYYWI REDGDYHNWT GCGNTLNLSH PAVVDYASAC LRYWVETCHV DGFRFDLAAV MGRTPEFRQD ALLFTAIQNC PVLSQVKLIA EPWDIAPGGY QVGNFPPLFA EWNDHFRDAA RRFWLHYDLP LGAFAGRFAA SSDVFKRNGR LPSAAINLVT AHDGFTLRDC VCFNHKXNEA NGEENRDGTN NNYSNNHGKE GLGGTLDLVE RRRDSIHALL TLLLSQGTP MLLAGDEHGH SQRGNNNAYC QDNQLTWLDW SQASSGLTAF TAALIHRLKR IPALMENRWW EEGDGNVRWL NRYAQLPLSTD EWQNGPKQLQ ILLSDRFLIA INATLEVTEI VLPAGEWHAI PPFAGEDNPV ITAVWQGPAH GLCVFQR
<b>Source</b>	E.coli
<b>Target Names</b>	glgX
<b>Protein Names</b>	Recommended name: Glycogen debranching enzyme EC= 3.2.1.- Alternative name(s): Glycogen operon protein GlgX
<b>Expression Region</b>	1-657
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.