



Recombinant Danio rerio Cyclin-dependent kinase 15 (cdk15)

Product Code	CSB-EP630193DIL
Abbreviation	cdk15
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q1RLU9
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MQNLRHAASE AFQRLGLKQR HLGYEELGEL DGVEKPQPHW FHTLQVRRRLR VQRGRSNDP MGGKSFQQEF QWKTGLQFGN ATSYLNLEKL GEGTYATVYK GISRINGHLV ALKVIHMKTE EGIPFTAIRE ASLLKGLKHA NIVLLHDIH TRESLTFVFE YVQTDLAQYM IQHPGGLHSY NIRLFMFQLL RGLSYIHGRR ILHRDLKPQN LLISYLGELK LADFGLARSK SIPCQTYSAE VVTLWYRPPD VLMGSTDYST ALDIWGAGCI FIEMLQGSPA FPGVADVFEQ LLKIWTVIGV PTEEIWPGVS DLPNYKPEWF LPCKPQQFRD VWKRLSQLPY KTEDLAQQML MMNPKDRISA QDALLHPYFN TLPPPLMHLR DTVSIFKVPG VRLESEARDI FSPSRRTKTP LAPLAKCW
Source	E.coli
Target Names	cdk15
Protein Names	Recommended name: Cyclin-dependent kinase 15 EC= 2.7.11.22 Alternative name(s): Cell division protein kinase 15
Expression Region	1-418
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.