



Recombinant Guanine nucleotide-binding protein alpha-5 subunit (gpa-5)

Product Code	CSB-EP630296CXY
Abbreviation	gpa-5
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q20701
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	GLVTCKQER EAELQNRQID TQIRIENQAN KRKIKMLLLG VTDSGKSTIV KQMRVNYLDG FNETEVVNAI FVIRNNIIDA FKNICNIILH SDITVTQEEK VLVKLFAYES GKIELMQEVD ELNVINAVSG YECIKQFFER FAFHPMVPDH IHYFYPNLDR IASSNYVPTA EDLIHMRQTT LGVHEISFDY TKHIIRLIDV GGQKTERRKW IHFFEGVTAV MFVCSLASFN QTTEEEP KAF VWESSLNKVQ NKVLVRSAGK AKVEKPGLIN RLDESVDL FK SIRENSFLKM SNFMLFLNKK DLLTKKLT KV VFSDYFPDYK KWITNDNSDV SVAEFIENMF REGLEPEKRM YAHLTQATVT ANIEGTFALC CDVIFGKNYE DTNLE
Source	E.coli
Target Names	gpa-5
Protein Names	Recommended name: Guanine nucleotide-binding protein alpha-5 subunit
Expression Region	2-385
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.