



Recombinant Homeobox protein ceh-39 (ceh-39)

Product Code	CSB-BP630477CXY
Abbreviation	ceh-39
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q22812
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	<p>MDFSNTYRNY GEVVDVPEDF LSDYVPTVKS EPEPAAIAPS VCEPPISTTR NYGTDFGRAS SMCGSSSSSS SSSYSSGSSP NYLTENGGDT EEP RFEELVS RKMSFDSENV AQVSNQMDYA RDTYADSRIP FASIDSN DAL LPYPVNNYVA ALPVPVPAKK SRGRTAPRTP KSLETSKPGR GVKRPASKEI DAKDLGIGDL SAHMNRQIGD DEELDTQLIA HRILDELKEQ CIPQASLAVK VLGRSQGTLS DLLRKPYPWG IMKNGRGTFQ RMANWLDLDP VVRRALCFMK KEDVARITGM AEPTPAKRAR QTTPSDERIR RFTFTQTQLD SLHTVFQQQD RPNREMQQAL SATLKLNRST VGNFFMNARR RLPKAAVPSQ NLVAEHAI DH NQPGPSHHY</p>
Source	Baculovirus
Target Names	ceh-39
Protein Names	Recommended name: Homeobox protein ceh-39
Expression Region	1-399
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>