



# Recombinant Homeobox protein ceh-32 (ceh-32)

<b>Product Code</b>	CSB-EP630482CXY-B
<b>Abbreviation</b>	ceh-32
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q23175
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Caenorhabditis elegans
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MFTPEQFTKV MSQLGNFSQL GQMFQPGNVA MLQALQANGA SSTPSLFPAM PSVIPSLAAP SSPTTSLNLT DQIVKTCEQL ETDGDVDGLF RFMCTIPPQK TQEVAGNEAF LRARALVCFH ASHFRELYAI LENNKFSKY HPKLQEMWHE AHYREQEKNR GKSLCAVDKY RVRKKYMPR TIWDGEQKTH CFKERTRSL RYWLKDPYP NPPKKKELAN ATGLTQMQVG NWFKNRRQRD RAAAKNKQN IIGVELKKTS SDMSDSDDDF EDSMTDSPSP IDEPKDLKSKS HIPKLSPTLL PKMATPFDMF AAAANPLMML NLNPALYMQF HNFFNTMRNP QIDEEENSET TVEVEADIEP PKKRKLSID EILNIKSEVS PSQCSPCSNE SLSPKRAVKT EEVKKEDDEA AEEDSRSVKS ETSSEDPKHSS PKSTTSQSE
<b>Source</b>	E.coli
<b>Target Names</b>	ceh-32
<b>Protein Names</b>	Recommended name: Homeobox protein ceh-32
<b>Expression Region</b>	1-439
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.