



Recombinant Sel1 repeat-containing protein 1 homolog (ZK20.4)

Product Code	CSB-MP630487CXY
Abbreviation	ZK20.4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q23450
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MSGKISEEAM QEQAERA EYV KNIGIEYRFG CYEEKRPESC QLLGEYMEAL QQNFQAAVAI FKKNC DENS Y PKSCFKY GMY MLAGKVGDDA SLSKMIRPMK IACDANIPQG CRYLSLVHWN GEKDRKADSE KAERYMRRAC ELEDGEACWL LSTWYMGNKE KFKTNAAGEA KELDRAELGS LSRDMDKSLQ FAIRACEMDI PQSCANVSRM FKLGDGVPKD LDQAKKYADR AKEIMETLRS NDTNAGFTG
Source	Mammalian cell
Target Names	ZK20.4
Protein Names	Recommended name: Sel1 repeat-containing protein 1 homolog Alternative name(s): Beta-lactamase hcp-like protein ZK20.4
Expression Region	1-249
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.