



# Recombinant Drosophila melanogaster DNA primase small subunit (DNApol-alpha50)

|                          |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
|--------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-EP630508DLU-B                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Abbreviation</b>      | DNApol-alpha50                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                                                                                                                                      |
| <b>Uniprot No.</b>       | Q24317                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Product Type</b>      | Recombinant Protein                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Immunogen Species</b> | Drosophila melanogaster (Fruit fly)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Purity</b>            | ≥85% (SDS-PAGE)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Sequence</b>          | MPEQVTDQEN QAPQQQTTAV HAYNPEVLQD MLPVYYRRLF PHLPFYRWLS<br>YGSSDAIFS NREISFTLQD DIYIRYLCFD TQAELEKEIC SRNPIKIDIG<br>PVMHSPKPNH RSIPGGLTPV QRELVFDIDM TDYDEVRTCC SGAGVCLKCW<br>KFMVLAARVL DVALREDFGF EHIIWIFSGR RGIHCWVCDY QARHLDGRGR<br>YAVAELNII TYASFAGGNS PRCSMGDRPH HSLKRALKIV EPMFEEIVLE<br>DQNLFGTPKG VTKLLNMVHD DAARGELESY MQKNLEDGAH SRLVWESFIK<br>YANSMRTSTT SAWSRKLKNI VAEIQLGLLY PRLDINVTRG FNHLLKAPFC<br>IHPATGKVCV PFSVSAVAKF DPTTVPTITQ LLHEINAFDD KSKSYMEAPE<br>DKSRIKDHKK TSMFKGVVVF EEFLRKLERS QKSASLQF |
| <b>Source</b>            | E.coli                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Target Names</b>      | DNApol-alpha50                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| <b>Protein Names</b>     | Recommended name: DNA primase small subunit EC= 2.7.7.- Alternative name(s): DNA polymerase subunit A DNA primase 50 kDa subunit Short name= dPRI50                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Expression Region</b> | 1-438                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Protein Length</b>    | full length protein                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.                                                                                                             |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life                                                                                                                                                                                                                                                                                       |



of lyophilized form is 12 months at -20°C/-80°C.