



Recombinant Uncharacterized ORF3 protein (ORF3)

Product Code	CSB-EP630662NAAM
Abbreviation	ORF3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q27YG8
Product Type	Recombinant Protein
Immunogen Species	Nora virus
Purity	≥85% (SDS-PAGE)
Sequence	MALKEEIFDQ NTTLFAVLDE NEVTEIKSIQ SSVTAVKTQL DQQKLQLDGL AKVVDNNQAR NEEQFVNINT TLVEMSSEVD KLTITTSQQA KQINTLATTL NELDQTTKDS LDTLNTTTES LKKQVLFNTD EITVLKVDVA TVTQKQQDVE HSLVTMKDEI GELHVSVNAN ANSIEALRTR IAAL EIRDVG PWVLKDRIYK FVINKPNGTT RYTTIYFFAD VYYSTGVRAA PTNSGTATS I LTITSLTTSY SLANVPVLKG VPYRVNGYFA NGNNIEDITG STSVIYDSM
Source	E.coli
Target Names	ORF3
Protein Names	Recommended name: Uncharacterized ORF3 protein
Expression Region	1-289
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.