



Recombinant Mouse Heat shock 70 kDa protein 1-like (Hspa1l)

Product Code	CSB-MP010823MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P16627
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	<p>MAANKGMAIG IDLGTTYSCV GVFQHGKVEI IANDQGNRTT PSYVAFTDTE RLIGDAAKNQ VAMNPQNTVF DAKRLIGRKF NDPVVQSDMK LWPQVINEA GKPKVMVSYK GEKKAFFYPEE ISSMVLTKMK ETAEAFGLHN VTNNAVITVPA YFNDSQRQAT KDAGVIAGLN VLRIINEPTA AAIAYGLDKG SHGERHVLIF DLGGGTFDVS ILTIDDGIFE VKATAGDTHL GGEDFDNRLV SHFVEEFKRK HKKDISQNKR AVRRLRTACE RAKRTLSSST QANLEIDSLY EGIDFYTSIT RARFEELCAD LFRGTLEPVE KSLRDAKMDK AKIHDIVLVG GSTRIPKVQK LLQDYFNDRD LNKSINPDEA VAYGAAVQAA ILMGDKSEKV QDLLLLDVAP LSLGLETAGG VMTVLIKRNS TIPTKQTQIF TTYSDNQPGV LIQVYEGERA MTRDNLLGR FDLTGIPPAP RGV PQIEVTF DIDANGILNV TAMDKSTGKA NKITITNDKG RLSKEEIERM VQEAERYKAE DEGQREKIAA KNALESYAFN MKSAVGDEGL KDKISESDKK KILDKCNEVL SWLEANQLAE KDEFDHHRKE LENMCPNPIIT KLYQSGCTGP TCTPGYTPGR AATGPTIEEV D</p>
Source	Mammalian cell
Target Names	Hspa1l
Protein Names	Recommended name: Heat shock 70 kDa protein 1-like Short name= Heat shock 70 kDa protein 1L Alternative name(s): Heat shock 70 kDa-like protein 1 Spermatid-specific heat shock protein 70
Expression Region	1-641
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Target Details	This gene encodes a 70kDa heat shock protein. In conjunction with other heat shock proteins, this protein stabilizes existing proteins against aggregation and mediates the folding of newly translated proteins in the cytosol and in organelles. The gene is located in the major histocompatibility complex class III region, in a cluster with two closely related genes which also encode isoforms of the 70kDa heat shock protein.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.