



Recombinant Escherichia coli Aliphatic sulfonates import ATP-binding protein SsuB (ssuB)

Product Code	CSB-YP631152EGW
Abbreviation	ssuB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q1RDS4
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain UTI89 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MNTARLNQGT PLLLNAVSKH YAENIVLNQL DLHIPAGQFV AVVGRSGGGK STLLRLLAGL ETPTAGDVLA GTTPLAEIQD DTRMMFQDAR LLPWKSVIDN VGLGLKGQWR DAARQALAAV GLENRAGEWP AALSGGQKQR VALARALIHR PGLLLLDEPL GALDALTRLE MQDLIVSLWQ EHGFTVLLVT HDVSEAVAMA DRVLLIEEGK IGLDLTVDIP RPRRLGSVRL AELEAEVLQR VMQRGHSEQP IRRHG
Source	Yeast
Target Names	ssuB
Protein Names	Recommended name: Aliphatic sulfonates import ATP-binding protein SsuB EC= 3.6.3.-
Expression Region	1-255
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.