



Recombinant Escherichia coli UDP-2,3-diacetylglucosamine hydrolase (IpxH)

| | |
|--------------------------|---|
| Product Code | CSB-EP631160EGW |
| Abbreviation | IpxH |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q1RF10 |
| Product Type | Recombinant Protein |
| Immunogen Species | Escherichia coli (strain UTI89 / UPEC) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MATLFIADLH LCVEEPAITA GFLRFLAGEA RKADALYILG DLFEAWIGDD DPNPLHHQMA AAIKAVSDSG VPCYFIHG NR DFLLGKRFAR ESGMTLLPEE KVLELYGRRV LIMHGDTLCT DDAGYQAFRA KVHKPWLQTL FLALPLFVRK RIAARMRANS KEANSSKSLA IMDVNQNAV V SAMEKHQVQW LIHGHTHRPA VHELIANQQP AFRVVLGAWH TEGSMVKVTA DDVELIHFPF |
| Source | E.coli |
| Target Names | IpxH |
| Protein Names | Recommended name: UDP-2,3-diacetylglucosamine hydrolase EC= 3.6.1.54 Alternative name(s): UDP-2,3-diacetylglucosamine diphosphatase |
| Expression Region | 1-240 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |