



Recombinant ADAM family mig-17 (mig-17)

Product Code	CSB-BP631334CXY
Abbreviation	mig-17
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q20930
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	≥85% (SDS-PAGE)
Sequence	ARREKQQSND ISFVKRKVQD GLKFSRVIKY TNETIQGMKT NFNSNKTQEL SLDVLVVADF LSYQAFLEMS NGDSHRAIHN LKEYLHALFE QTKIYDGIS FGNETLHMVF AGTWIATQER DCPLWISWAE EEEERVLNEE IRRLEEKERD LNSTFVDDTF FMNSTDSDNS STDALISSDM PKKLRKFVDI TLEEMQENNS TEMTLKIDSK KAIDKFTIWL KEQTGLPRHE HAVLITKFDL ISINGNSATQ GMAYVGNICE NGDSSSVVED IGAGLTSLIM AHEIGHSLGA LHDGAYETAE CDSNDNYLMA VAVSGSADRQ SFLNSRRMSN CSINSIIENL KEPTANCVKK WTKKKGKDVS QKDFIKKPGE LVKITRQCQV AFGPTFIPCL HIGYFHEQSI CERIWCSDGE SDECQTLNYF PAFDGTECGY NMWCLEGSCV QNTKKWMDCK DINSKTCSKY STSKLKHYCK SKDFREICCR TCAKKGKIY
Source	Baculovirus
Target Names	mig-17
Protein Names	Recommended name: ADAM family mig-17 EC= 3.4.24.- Alternative name(s): A disintegrin and metalloproteinase mig-17 Abnormal cell migration protein 17
Expression Region	21-509
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.