



# Recombinant ADAM family mig-17 (mig-17)

<b>Product Code</b>	CSB-EP631334CXY
<b>Abbreviation</b>	mig-17
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q20930
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Caenorhabditis elegans
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	ARREKQQSND ISFVKRKVQD GLKFSRVIKY TNETIQGMKT NFNSNKTQEL SLDVLVVADF LSYQAFLEMS NGDSHRAIHN LKEYLHALFE QTKIYDGIS FGNETLHMVF AGTWIATQER DCPLWISWAE EEEERVLNEE IRRLEEKERD LNSTFVDDTF FMNSTDSDNS STDALISSDM PKKLRKFVDI TLEEMQENNS TEMTLKIDSK KAIDKFTIWL KEQTGLPRHE HAVLITKFDL ISINGNSATQ GMAYVGNICE NGDSSSVVED IGAGLTSLIM AHEIGHSLGA LHDGAYETAE CDSNDNYLMA VAVSGSADRQ SFLNSRRMSN CSINSIIENL KEPTANCVKK WKTKKGKDVS QKDFIKKPGE LVKITRQCQV AFGPTFIPCL HIGYFHEQSI CERIWCSDGE SDECQTLNYF PAFDGTECGY NMWCLEGSCV QNTKKWMDCK DINSKTCSKY STSKLKHYCK SKDFREICCR TCAKKGKIY
<b>Source</b>	E.coli
<b>Target Names</b>	mig-17
<b>Protein Names</b>	Recommended name: ADAM family mig-17 EC= 3.4.24.- Alternative name(s): A disintegrin and metalloproteinase mig-17 Abnormal cell migration protein 17
<b>Expression Region</b>	21-509
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.