



Recombinant Tropomyosin isoforms a/b/d/f (lev-11)

Product Code	CSB-MP631512CXY
Abbreviation	lev-11
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q22866
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	≥85% (SDS-PAGE)
Sequence	MDAIKKKMQA MKIEKDNALD RADAAEEKVR QITEKLERVE EELRDTQKKM TQTGDDLDKA QEDLSAATSK LEEKEKTVQE AEAEVASLNR RMTLLEEELE RAEERLKIAT EKLEEATHNV DESERVRKVM ENRSLQDEER ANTVEAQLKE AQLLAEEADR KYDEVARKLA MVEADLERAE ERAEAGENKI VELEELRVV GNNLKSLEVS EEKALQREDS YEEQIRTVSS RLKEAETRAE FAERSVQKLQ KEVDRLEDEL VHEKERYKTI SEELDSTFQE LSGY
Source	Mammalian cell
Target Names	lev-11
Protein Names	Recommended name: Tropomyosin isoforms a/b/d/f Alternative name(s): Levamisole resistant protein 11
Expression Region	1-284
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.