



Recombinant *Drosophila melanogaster* Putative odorant-binding protein A10 (a10)

Product Code	CSB-BP631664DLU
Abbreviation	a10
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q27377
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	L PHPPATSPSP MMERMVEQAY DDKFDNVDLD EILNQERLLI NYIKCLEGTG PCTPDAKMLK EILPDAIQTD CTKCTEKQRY GAEKVTRHLI DNRPTDWERL EKIYDPEGTY RIKYQEMKSK ANEEP
Source	Baculovirus
Target Names	a10
Protein Names	Recommended name: Putative odorant-binding protein A10 Alternative name(s): Antennal protein 10 OS-D protein
Expression Region	30-155
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.