



Recombinant Danio rerio Poly (U)-specific endoribonuclease-C (endouc)

Product Code	CSB-BP631873DIL
Abbreviation	endouc
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q1LUM3
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	SSQTVNQE LSNIFNELWK LDVNRMEPLT NYNISLQGKA GYIPQGSTNV VDHASSPLFV NVDEAKLSSI TTYARFMKLL DNYERSTGVA ERVTAAEEVTE NNSFLDAILE TAVMKRAHQY LIGKGKSRSD LRQFKSQLYY MWFRLYHRER NGGEDSSGFE HVFVGETKFG REIMGLHNWV QFYLQEKQNL LDYKGYKARA NDVPDADDHV LNVQFSWHGL VKPVASAFVG VSPEFEMAVF TILFLTSTEK TTTAVVNLDE YQLEMVVHRH GRCIGTAYPK LLSSNNRHM
Source	Baculovirus
Target Names	endouc
Protein Names	Recommended name: Poly(U)-specific endoribonuclease-C EC= 3.1.-.- Alternative name(s): Protein endoU-C Uridylate-specific endoribonuclease-C
Expression Region	23-309
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.