



Recombinant Danio rerio UPF0524 protein C3orf70 homolog A (si:dkey-24b15.2)

Product Code	CSB-EP631880DIL-B
Abbreviation	si:dkey-24b15.2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q1LY84
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MAASGGQKSE KLDEAQALAR SCAGRPDFLP CDGLSICATH SHGKCFKLHW CCHLGWCHCK YVYQPMTSVC QLPSTAVPVA PSGHTHTMNL SISLAERFLR TAPKFQAPPK PESPKFCVIS DLFVDDYMKV RINGKMCYVQ RPPPPTPNPT HQPQTAAPQP VPQRSKNSHQ VGPTVKQGQA KEKISAPKMD HCSSPSSSED SGINALGLHY MESCEDESDCV DDDDEEEEDD ELSTDGNSSP GSFWDQDECT LLSPSKSIVE IIEKIETTV
Source	E.coli
Target Names	si:dkey-24b15.2
Protein Names	Recommended name: UPF0524 protein C3orf70 homolog A
Expression Region	1-269
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.