



Recombinant Escherichia coli UPF0207 protein yfbR (yfbR)

Product Code	CSB-YP632158EGW
Abbreviation	yfbR
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q1R9C4
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain UTI89 / UPEC)
Purity	≥85% (SDS-PAGE)
Sequence	MKQSHFFAHL SRLKLINRWP LMRNVRTENV SEHSLQVAMV AHALAAIKNR KFGGNVNAER IALLAMYHDA SEVLGTDLPT PVKYFNSQIA QEYKAIEKIA QQKLVDMVPE ELQDIFAPLI DEHAYSDEEK SLVKQADALC AYLKCLEELA AGNNEFLLAK TRLEATLEAR RSQEMDYFME VVPSFHLSL DEISQDSPL
Source	Yeast
Target Names	yfbR
Protein Names	Recommended name: UPF0207 protein yfbR
Expression Region	1-199
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.