



# Recombinant Human Zinc finger protein 23 (ZNF23)

<b>Product Code</b>	CSB-BP026597HU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P17027
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MLENYGNVAS LGFPLLKPAV ISQLEGGSEL GGSSPLAAGT GLQGLQTDIQ TDNDLTKEMY EGKENVSFEL QRDFSQETDF SEASLLEKQQ EVHSAGNIKK EKSNTIDGTV KDETSPVEEC FFSQSSNSYQ CHTITGEQPS GCTGLGKSIS FDTKLVKHEI INSEERPFKC EELVEPFRCD SQLIQHQENN TEEKPYQCSE CGKAFSINEK LIWHQRLHSG EKPFKVECG KSFSYSSHYI THQTIHSGEK PYQCKMCGKA FSVNGSLSRH QRIHTGEKPY QCKECCNGGFS CSSAYITHQR VHTGEKPYEC NDCGKAFNVN AKLIQHQRH TGEKPYECNE CGKGFRCSSQ LRQHQSHTG EKPYQCKECCG KGFNNNTKLI QHQRIHTGEK PYECTECGKA FSVKGLIQH QRIHTGEKPY ECNECGKAFR CNSQFRQHLR IHTGEKPYEC NECGKAFSVN GKLMRHQRH TGEKPFECNE CGRCFTSKRN LLDHHRHTG EKPYQCKECCG KAFSINAKLT RHQRHTGEK PFKCMECEKA FSCSSNYIVH QRIHTGEKPF QCKECCGKAFH VNAHLIRHQR SHTGEKPFRC VECGKGFSS SDYIIHQTVH TWKKPYMCSV CGKAFRFSFQ LSQHQSVMHSE GKS
<b>Source</b>	Baculovirus
<b>Target Names</b>	ZNF23
<b>Protein Names</b>	Recommended name: Zinc finger protein 23 Alternative name(s): Zinc finger protein 359 Zinc finger protein 612 Zinc finger protein KOX16
<b>Expression Region</b>	1-643
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.