



Recombinant Human Zinc finger protein 23 (ZNF23)

Product Code	CSB-MP026597HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P17027
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MLENYGNVAS LGFPLLKPAV ISQLEGGSEL GGSSPLAAGT GLQGLQTDIQ TDNDLTKEMY EGKENVSFEL QRDFSQETDF SEASLLEKQQ EVHSAGNIKK EKSNTIDGTV KDETSPVEEC FFSQSSNSYQ CHTITGEQPS GCTGLGKSIS FDTKLVKHEI INSEERPFKC EELVEPFRCD SQLIQHQENN TEEKPYQCSE CGKAFSINEK LIWHQRLHSG EKPFKVECG KSFSYSSHYI THQTIHSGEK PYQCKMCGKA FSVNGSLSRH QRIHTGEKPY QCKECCNGGFS CSSAYITHQR VHTGEKPYEC NDCGKAFNVN AKLIQHQRH TGEKPYECNE CGKGFRCSSQ LRQHQSHTG EKPYQCKECCG KGFNNNTKLI QHQRIHTGEK PYECTECGKA FSVKGLIQH QRIHTGEKPY ECNECGKAFR CNSQFRQHLR IHTGEKPYEC NECGKAFSVN GKLMRHQRH TGEKPFECNE CGRCFTSKRN LLDHHRHTG EKPYQCKECCG KAFSINAKLT RHQRHTGEK PFKCMECEKA FSCSSNYIVH QRIHTGEKPF QCKECCGKAFH VNAHLRHQR SHTGEKPFRC VECGKGFSS SDYIIHQTVH TWKKPYMCSV CGKAFRFSFQ LSQHQSVMHSE GKS
Source	Mammalian cell
Target Names	ZNF23
Protein Names	Recommended name: Zinc finger protein 23 Alternative name(s): Zinc finger protein 359 Zinc finger protein 612 Zinc finger protein KOX16
Expression Region	1-643
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.