



Recombinant Elongation factor Ts, mitochondrial (tsfm-1)

Product Code	CSB-BP632371CXY
Abbreviation	tsfm-1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q20819
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	AE AEKKVSKEAL MALRKKTGYS YVNCRKALIQ FGENDMDSAV KWLKEAAAKE GWAKAAKLGTV RVTNGLVSV VTDNSTAAVV ELSCETDFVA RSGAFKDLLS NISNSVLAKA KPQSISSGSK LQEFTYDLGD LTDSGKNMR EVLSSLSIGKL GENMTVRRVK AFKAPEGTTL FGASHPKDGT DDIPMGRFIS LIALNQSSPG SISSQQLAGQ ICQHIIGMSP ESLGEEAESV KTQEGLRSQE GHDPNADPVV VTNIDESETA LLRQAFMLNP SQSVHEYLKS HNANILDFVR VELGSE
Source	Baculovirus
Target Names	tsfm-1
Protein Names	Recommended name: Elongation factor Ts, mitochondrial Short name= EF-Ts Short name= EF-TsMt
Expression Region	19-316
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.