



Recombinant Probable 26S proteasome regulatory subunit rpn-6.1 (rpn-6.1)

Product Code	CSB-EP632373CXY
Abbreviation	rpn-6.1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q20938
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MRETSSREDT NNIGKAPEMS GGTIMDTMTS LPHQNDQNV I RHLTNLVKSP ASGDDDIKKK EDSIMELGNI LAQNKQTEEL RNMIEQTRPF LVSLGKAKAA KLVRDLVDLC LKIDDQDGI KVGLVKECIQ WATEQNRFL RQTLTARLVR LYNDLQRYTQ ALPLAADLIR ELKKVDDKDV LVEVELEESK AYYNLSNIGR ARASLTGART TANAIYVNPR MQAALDLQSG ILHAADEKDF KTAFSYFYEA FEGYDSVDEK VSALTALKYM LLCKVMLDLP DEVNSLLSAK LALKYNGSDL DAMKAIAAAA QKRSLKDFQV AFGSFPQELQ MDPVVRKHFH SLSERMLEKD LCRIIEPYSF VQIEHVAQQI GIDRSKVEKK LSQMILDQKL SGSLDQGEEM LIVFEIAVPD EAYQTALDTI HAMGEVVDAL YSNASKIN
Source	E.coli
Target Names	rpn-6.1
Protein Names	Recommended name: Probable 26S proteasome regulatory subunit rpn-6.1
Expression Region	1-438
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.