



Recombinant Chondroitin proteoglycan 3 (cpg-3)

Product Code	CSB-BP632426CXY
Abbreviation	cpg-3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q21771
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	≥85% (SDS-PAGE)
Sequence	HPA DPIRAKRDVS ASEDEFSGDS SGEISGESSG EASGEASGEA SGEASGEASG ESSGETSGES SGDEETSGEG SGEEGSGDTS PVVPVDELTL QQLETLNTYA QQVQAESQKL IHQANFVITE MTALSANAQN LGILSNIVLA NSQMVLD SAR LSLNETETET GTSAPATCVS SAVCYGDSGC GSGKCIGALA GTCNCNSCVF GWPCQEDSAC GGFNGACNSI TATCDCFAAY TKNLTLAEA LTSFCNVETC NGAEDNVEKC HGLPCNYGFC VC
Source	Baculovirus
Target Names	cpg-3
Protein Names	Recommended name: Chondroitin proteoglycan 3
Expression Region	18-292
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.