



Recombinant Nuclear hormone receptor family member nhr-11 (nhr-11)

Product Code	CSB-MP632557CXY
Abbreviation	nhr-11
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q23294
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MGPLCAVCES PTAFTLHFGG RCCKACAAFF RRTIALDLKY ECAADDPCEI HFKFPGMRLV CRECRLRKCY SAGMRSELVR SKRENFACR RKDSRNNSDA APNSNSPSTR QSSSPEEMDD WSFQMFEKPK KIEDLPLTPS ISHPPLPTQL MPEESSRTSS FDGGYCSSYP SSSAATHPSP PGMLYSIENN SILQYYHSME TGLCSKRRIM YTNMGMDFIL DTHANLQCPF TVNDLRPHDY RNFRGMLRHD FVMLFDYATR FPEFNSFTSH EKNMFYRQIV AVDFILSSAY YTAKLGQAHR QMVLNNGEYL NMDPLMSGN EIDARRYFDC DEDFSKYRAL MPMHIAIWE SIVPFSKLVN TFEFCLLKA LTVWQATYFK LTENGREKCR RQRNIIIGFL SKMCHSPGGG GEHRVGELLM SMNYLRESAQ KLTTSYVMLT VFNVLNCDSM LHEMLNFQY
Source	Mammalian cell
Target Names	nhr-11
Protein Names	Recommended name: Nuclear hormone receptor family member nhr-11
Expression Region	1-459
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.