



Recombinant Drosophila melanogaster Cell death protein W (W)

Product Code	CSB-EP632574DLU
Abbreviation	W
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q24106
Product Type	Recombinant Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MAVPFYLP EG GADDVASSSS GASGNSSPHN HPLPSSASSS VSSSGVSSAS ASSASSSSSA SSDGASSAAS QSPNTTTSSA TQTPMQSPLP TDQVLYALYE WVRMYQSQQS APQIFQYPPP SPSCNFTGGD VFFPHGHPNP NSNPHPRTPR TSVSFSSGEE YNFFRQQQPQ PPHSYAPST P QPMPPQSAP PMHCSHSYPQ QSAHMMPHHS APFGMGGTY Y AGYTPPPTPN TASAGTSSSS AAFGWHGHPH SPFTSTSTPL SAPVAPKMRL QRSQSDAARR KRLTSTGEDE REYQSDHEAT WDEFGDRYDN FTAGRERLQE FNGRIPPRKK KSSNSHSSSS NNPVCHTDSQ PGGTQAESG AIHGHSQQR QVERERQKAK AEKKK PQSFT WPTVVTVFVL AMGCGFFAAR
Source	E.coli
Target Names	hid
Protein Names	Recommended name: Cell death protein W Alternative name(s): Protein head involution defective Protein wrinkled
Expression Region	1-410
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.