



Recombinant Callithrix jacchus Interleukin-3 (IL3)

Product Code	CSB-EP632965CYL-B
Abbreviation	IL3
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q28334
Product Type	Recombinant Protein
Immunogen Species	Callithrix jacchus (White-tufted-ear marmoset)
Purity	>85% (SDS-PAGE)
Sequence	AA PTQTMLP LTT QVNC SNLREE IVTLLNQ PPL ASSNFNNLNR EDQRILMKPN LRKPNLEAFQ KAVKSLQNA A AIESNLKDL P VCLPTATNAA TQHPIRIKDG DWNDFQM K LK FYLKTLEIKQ PQ
Source	E.coli
Target Names	IL3
Protein Names	Recommended name: Interleukin-3 Short name= IL-3 Alternative name(s): Hematopoietic growth factor Mast cell growth factor Short name= MCGF Multipotential colony-stimulating factor P-cell-stimulating factor
Expression Region	19-142
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>