



# Recombinant Danio rerio LIM/homeobox protein Lhx9 (Ihx9)

<b>Product Code</b>	CSB-MP633128DIL
<b>Abbreviation</b>	Ihx9
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q1LWV4
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Danio rerio (Zebrafish) (Brachydanio rerio)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MEVVGCKAKE SSCTLRPAAM LFHGISGDHI QGIMEEMERR SKSESRLAKT VQMNGRETTM PSMSPEKPAL CAGCGGKISD RYYLLAVDKQ WHLRCLKCCE CKLALESELT CFAKDGSIIYC KEDYYRRFSV QRCARCHLGI SASEMVMRAR DSVYHLSCFT CTTCNKLTLLT GDHFGMKDNL VYCRVHFETL IQGEYHPQLN YAELAAKGGG LALPYFNGTG TVQKGRPRKR KSPAMGIDIG SYSSGCNEND ADHLDRDQQP YPPSQKTKRM RTSFKHHQLR TMKSYFAINH NPDAKDLKQL AQKTGLTKRV LQVWFQARA KFRRNVLRQE NGGVDKADGT SLPPPSSDSG ALTPPSTATT LTDLTNPSIT VVTSVTSSLD SHESGSPQQT TLTNLF
<b>Source</b>	Mammalian cell
<b>Target Names</b>	Ihx9
<b>Protein Names</b>	Recommended name: LIM/homeobox protein Lhx9 Short name= LIM homeobox protein 9
<b>Expression Region</b>	1-396
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.