



Recombinant Oncorhynchus mykiss Tumor necrosis factor receptor type 1-associated DEATH domain protein (tradd)

Product Code	CSB-BP633137OEI
Abbreviation	tradd
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q1M161
Product Type	Recombinant Protein
Immunogen Species	Oncorhynchus mykiss (Rainbow trout) (Salmo gairdneri)
Purity	>85% (SDS-PAGE)
Sequence	MVALSVEVGP WTGCAVLFLR SCPAVNLLSL YKDQQGKFSV FKVLKLTLD SAGGLEGYEI LKLHDADPLL GVEVKFVDVA ACRRFLQSYG SGAVQQSLSQ HTCROLLHIPQ QELALETQLK AGTHTLDFCL DDLELCLQHI HQSQPERLRD DEIAELDQQL QSQUALGHIPQ PPTLTQEEPP VPSNCFRLFQK RVFEDRMLAA GDLQRFSNGV GRDWRKVGRA LGKNCRALKG PAIDNLAYEY EREGLYEQAY QLLSRFIQAE GRAARLGRLV RALEDSKLTLS LAENILDIQP RD
Source	Baculovirus
Target Names	tradd
Protein Names	Recommended name: Tumor necrosis factor receptor type 1-associated DEATH domain protein Short name= TNFR1-associated DEATH domain protein Alternative name(s): TNFRSF1A-associated via death domain protein
Expression Region	1-292
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.