



Recombinant Escherichia coli Uncharacterized protein YjeA (yjeA)

Product Code	CSB-BP633363EGW
Abbreviation	yjeA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q1R3A2
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli (strain UTI89 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MSETASWQPS ASIPNLLKRA AIMAEIRRF ADRGVLEVET PCMSQATVTD IHLVPFETRF VGPGHSQGMN LWLMTSPEYH MKRLLVAGCG PVFQLCRSFR NEEMGRYHNP EFTMLEWYRP HYDMYRLMNE VDDLLQQVLD CPAAESLSYQ QAFLRYLEID PLSADKTQLR EVAAKLDLSN VADTEEDRDT LLQLLFTFGV EPNIGKEKPT FVYHFPASQA SLAQISTEDH RVAERFEVYY KGIELANGFH ELTDAREQQQ RFEQDNRKRA ARGLPQHPID QNLIEALKVG MPDCSGVALG VDRLVMLALG AETLAEVIAF SVDRA
Source	Baculovirus
Target Names	epmA
Protein Names	Recommended name: Uncharacterized protein YjeA
Expression Region	1-325
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.