



# Recombinant Escherichia coli Isocitrate dehydrogenase kinase/phosphatase (aceK)

<b>Product Code</b>	CSB-YP633366EGW
<b>Abbreviation</b>	aceK
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q1R3T1
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli (strain UTI89 / UPEC)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MPRGLELLIA QTILQGFDAQ YGRFLEVTSG AQQRFEQADW HAVQQAMKNR IHLYDHHVGL VVEQLRCITN GQSTDAAFLL RVKEHYTRLL PDYPRFEIAE SFFNSVYCRL FDHRSITPER LFIFSSQPER RFRTIPRPLA KDFHPDHGWE SLLMRVISDL PLRLRWQNKS RDIHYIVRHL TETLGTDNLA ESHLQVANEL FYRNKAAWL V GKLITPSGTL PFLLIHQTD DGELFIDTCL TTTAEASIVF GFARSYFMVY APLPAALVEW LREILPGKTT AELYMAIGCQ KHAKTESYRE YL VYLQGCNE QFIEAPGIRG MVMLVFTLPG FDRVFKVIKD RFAPQKEMSA AHVRACYQLV KEHDRVGRMA DTQEFENFVL EKRHISPALM ELLLQEAEEK ITDLGEQIVI RHLYIERRMV PLNIWLEQVE GQQLRDAIEE YGN AIRQLAA ANIFPGDMLF KNFGVTRHGR VVFYDYDEIC YMTEVNFVDI PLPRYPEDEL ASEPWYSVSP GDVFPEEFRH WLCADPRIGP LFEEMHADLF RADYWRALQN RIREGHVEDV YAYRRRQRFS VRFV
<b>Source</b>	Yeast
<b>Target Names</b>	aceK
<b>Protein Names</b>	Recommended name: Isocitrate dehydrogenase kinase/phosphatase Short name= IDH kinase/phosphatase Short name= IDHK/P EC= 2.7.11.5 EC= 3.1.3.-
<b>Expression Region</b>	1-574
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at  $-20^{\circ}\text{C}/-80^{\circ}\text{C}$ . The shelf life of lyophilized form is 12 months at  $-20^{\circ}\text{C}/-80^{\circ}\text{C}$ .