



Recombinant Putative carbonic anhydrase-like protein 1 (cah-1)

Product Code	CSB-BP633623CXY
Abbreviation	cah-1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q20781
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	≥85% (SDS-PAGE)
Sequence	QNYQW SYDSDVFGGP HFWGLVEKDW WMCKKGRLQS PIDIQPDRLL FDASVKPVSY AELNFFLIY LTISLKLMRY IFTCAHLFAC NKKLGTWNRI FGANIAKFRA WVVYIIIFHA QVPQLQVVSE FVNTGQMVRV RIGYSSKKPS VNITSGPLYG YRYRVQRIDF HMGRKNENGS EHTINGRRFP MEVQLVAYNT DLYPNFTSAS KSPHGIAILS VLVDGFPETN QELIKLIAT ASISYKDQRV QLADFEPWRL LPFTRDIITY EGSLTSPGCH ETVTWILNQ PIFIKKEHFE EWSHLYLSME GAEKVPVAPN FRKIQETNNR LVRTNIQHKV WIFLS
Source	Baculovirus
Target Names	cah-1
Protein Names	Recommended name: Putative carbonic anhydrase-like protein 1
Expression Region	26-365
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.