



Recombinant UPF0046 protein K07C11.7 (K07C11.7)

Product Code	CSB-MP633659CXY
Abbreviation	K07C11.7
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q21268
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	KSQRIDPD AENELWDSIK HTKVQNIVEK NAIGEINSPP TDGTPYLKVV CISDTHEQLH NTVPDGDVL IHAGDFTNNG KREELIKFNE EMTRFPHKYK LVVAGNHELG FDHDENQGER QDADKGLGTE DGYNILTNVT YLQDKGVTID GVTFFGSSYH PLRGFPFYRN RAEQLAECWK AVPNDTNVLI THTPPLGYLD QFGDERWGCR DLLKTVERIQ PAYHIFGHVH EQHGVLSNGN TTFINAAQCN KGNQIQTRPI VFYIPKKTLS
Source	Mammalian cell
Target Names	K07C11.7
Protein Names	Recommended name: UPF0046 protein K07C11.7
Expression Region	23-290
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.