



Recombinant Probable voltage-dependent anion-selective channel (R05G6.7)

Product Code	CSB-MP633679CXY
Abbreviation	R05G6.7
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q21752
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MAPPTFADLG KSAKDLFNKG YNFGFLKIDS TTRAGDNKEV EFKSAASHNI GSGKLGGNLD VKYKIPQYGI TLTEKWNTEN QLGTVIEVNE QFGRGLKVTL DSLYAPHAGK RSGKVKLDWA LPTARVTADV GVTSAPVINA AGVFSRDGWL IGAAATFDSS SNKLAATSLA FGHSTPQYTL HSFVINSTDF GASLYHKVAS NVEVGTQLGW KVGNGADYA LATKYAPSRD LTVRAKVNSS SQVAVAATHS LSPALKLTLS TQFNLAANDA HKFGLGLEFD PSN
Source	Mammalian cell
Target Names	vdac-1
Protein Names	Recommended name: Probable voltage-dependent anion-selective channel
Expression Region	1-283
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.