



Recombinant Kinetochore protein Nuf2 homolog (him-10)

Product Code	CSB-MP633684CXY
Abbreviation	him-10
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q21952
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MSNVVLIVYD PRMISKYLGQ KLHMGLVADD IIKPTAEIAQ QIFANFVRLV LNVSESSLTT LPLSANCDYD PELHKKSIPI IILFQCMKAF IKDNSGNKLD LTMCDLVTPA KHEHRFRKLT SFLVDFLKLH ELATPAFNEI SEEFSDRKFE MEKIREELLE AEKKKNDLLA KQSIRKRHEH ELINEQSNK AELKNVVNEY TETRQINEEL DKQKEEAILH IQALEKEMLT GKKTIEHLNE EVLTSPEQLK QEMEERKRHI EELRDCLESS KKGLQAKLEA REICINSEKN VPVIEKIHQ WTEVREVIID LIDVESENLR KLKEMEEQLD FMMKEMETAQ KRLVEQSETH EQLRIHTQK SEERQRRIEE ITEQIANLKT SQPDVSQEIA KKKQELLALK NAHSETISQI TNSCQDAVAK FAKLNAMFKE TQKVAFEKNT AAAREMERLK SSLTGRLSD YTFGSSTIDA GENTENCDPQ PNDSSFSVFK
Source	Mammalian cell
Target Names	him-10
Protein Names	Recommended name: Kinetochore protein Nuf2 homolog Short name= CeNuf2 Alternative name(s): Kinetochore protein him-10 Protein high incidence of males 10
Expression Region	1-490
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.