



Recombinant Probable DNA replication complex GINS protein PSF1 (R53.6)

Product Code	CSB-MP633791CXY
Abbreviation	R53.6
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q22019
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	≥85% (SDS-PAGE)
Sequence	MSSGDQNRGG VADKALQLVL EMKRNPDLVLP PYNTELVLRQC YQKIDELFQK NAAVVEKIRA GLPHDSTLLQ PRLAAMCHIR RCMMAYVNER KNRIRSFRWK YGGALPASVR NALCDAEIQF FNEYSSTLAR FQSNLGEQGV NLLLHSAPPK SLFVQVRALE DYGEFETSDG TQVQLSKDSL HSLPRQDCM LIRQGVLELV H
Source	Mammalian cell
Target Names	psf-1
Protein Names	Recommended name: Probable DNA replication complex GINS protein PSF1 Alternative name(s): GINS complex subunit 1
Expression Region	1-201
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.