



Recombinant Probable splicing factor, arginine/serine-rich 2 (rsp-2)

Product Code	CSB-BP633809CXY
Abbreviation	rsp-2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q23120
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	≥85% (SDS-PAGE)
Sequence	MVRVYIGRLP NRASDRDVEH FFRGYGKLSD VIMKNGFGFV DFQDQRDADD AVHDLNGKEL CGERVILEFP RRKVGYNNEER SGSGFRGREP TFRKGGERQF SNRYSRPCST RFRLVIDNLS TRYSWQDIKD HIRKLGIEPT YSEAHKRNVN QAIVCFTSHD DLRDAMNKLQ GEDLNGRKLK CTDETRDRSR SRSPRRRSRS RSPTSRSPPP ARRRSPGSDR SDRKSRSASP KKRSDKRARS ESKSRSRSGG RRSRSNSPPN RSPSPKKRRD NSSPRSGSAS P
Source	Baculovirus
Target Names	rsp-2
Protein Names	Recommended name: Probable splicing factor, arginine/serine-rich 2 Alternative name(s): CeSRp40 RNA-binding protein srp-4
Expression Region	1-281
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.