



# Recombinant Drosophila melanogaster Dual specificity mitogen-activated protein kinase kinase dSOR1 (Dsor1)

<b>Product Code</b>	CSB-YP634507DLU
<b>Abbreviation</b>	Dsor1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q24324
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Drosophila melanogaster (Fruit fly)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MSKNKLNVLV PPVNTEATVA AATVAPTPPF KTPSGTDTHS LLGKPKTSID ALTETLEGLD MGDTERKRIK MFLSQKEKIG ELSDEDLEKL GELGSGNGGV VMKVRHTHTH LIMARKLIHL EVKPAIKKQI LRELKVLHEC NFPHIVGFYD AFYSDGEISI CMEYMDGGSL DLILKRAGRI PESILGRITL AVLKGLSYLR DNHAIHRDV KPSNILVNSS GEIKICDFGV SGQLIDSMAN SFVGTRSYMS PERLQGTHYS VQSDIWSLGL SLVEMAIGMY PIPPPNTATL ESIFADNAEE SGQPTDEPRA MAIFELLDYI VNEPPKLEH KIFSTEFKDF VDICKKQPD ERADLKTLLS HPWIRKAELE EVDISGWVCK TMDLPPSTPK RNTSPN
<b>Source</b>	Yeast
<b>Target Names</b>	Dsor1
<b>Protein Names</b>	Recommended name: Dual specificity mitogen-activated protein kinase kinase dSOR1 Short name= Downstream of RAF Short name= MAPKK EC= 2.7.12.2
<b>Expression Region</b>	1-396
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.