



Recombinant Abnormal cell lineage protein 44 (lin-44)

Product Code	CSB-BP634642CXY
Abbreviation	lin-44
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q27886
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	TNEIP TISGAPAGKI VQPPKPNILK QGCPSDLLHS RALRSIQLAC RTHPATVISA FEGVQEGLQN CANRLRFQQW DCSEAGNIMH DPPLLRQGFR ESSLIWALSS ASAAWGVATA CAQGWIDDCA CNNQMGQNEY EFGGCTHGVQ HGITASRKLK TKVGAVNTLL RKVEKHNLKA GRLAIKKTLL SSCKCHGVSG SCQQKTCWKR TATLEHITDY LVEKYARAKL YTDDSVVKTT DLIYLEASPD VCKAKSVAGR VCAWRNETHT QGDCDRGCCG NGFSIRHEVV RVKCDCEFWW CCNLVCKDCI QHRWISTCNG TPPKSLIF
Source	Baculovirus
Target Names	lin-44
Protein Names	Recommended name: Abnormal cell lineage protein 44 Alternative name(s): Wnt protein
Expression Region	26-348
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.